



Complete System

The Medicina SUPERLOCK[®] system is a complete range of equipment which has been developed for the safe delivery of enteral formula, feeds, medicines and flushes.



SUPERLOCK[®]

Safe & cost effective

The SUPERLOCK[®] connectors are bigger than the standard Luers and it is mechanically impossible to connect any device in the SUPERLOCK[®] system to any intravenous device under any circumstances.




Syringes

Complete range of colour coded Enteral Safety Syringes.



Nasogastric Feeding Tubes

A comprehensive range of sizes of nasogastric feeding tubes in both polyurethane or silicone.



Accessories

A wide range of accessories for medicines and feeds.

Pump

Dedicated enteral pump



Enteral Syringe Pump

P900	Enteral Syringe Pump	1
------	----------------------	---

Syringes

MD01	1ml SUPERLOCK Syringe	100
MD25	2.5ml SUPERLOCK Syringe	100
MD05	5ml SUPERLOCK Syringe	100
MD10	10ml SUPERLOCK Syringe	100
MD20	20ml SUPERLOCK Syringe	80
MD30	30ml SUPERLOCK Syringe	80
MD60	60ml SUPERLOCK Syringe	60
MD60B	60ml SUPERLOCK Syringe Cath.	55

Enteral Extension Sets

MF150Y	150cm Extension Set Female	50
MF150MM	150cm Extension Set Male	50
MF075	75cm Extension Set	50
MF150	150cm Extension Set	50

Enteral Drainage Bags

MED50M	50ml Enteral Drainage Bag	50
MED100M	100ml Enteral Drainage Bag	50
MED250M	250ml Enteral Drainage Bag	50
MED500M	500ml Enteral Drainage Bag	50

Enteral Feeding Tubes

MG4/40	4FR x 40cm Polyurethane	50
MG5/40	5FR x 40cm Polyurethane	50
MG5/50	5FR x 50cm Polyurethane	50
MG5/80	5FR x 80cm Polyurethane	50
MG5/90	5FR x 90cm Polyurethane	50
MG6/40	6FR x 40cm Polyurethane	50
MG6/50	6FR x 50cm Polyurethane	50
MG6/60	6FR x 60cm Polyurethane	50
MG6/80	6FR x 80cm Polyurethane	50
MG8/40	8FR x 40cm Polyurethane	50
MG8/50	8FR x 50cm Polyurethane	50
MG8/80	8FR x 80cm Polyurethane	50
MG8/90	8FR x 90cm Polyurethane	50
MG10/80	10FR x 80cm Polyurethane	50

Silicone Enteral Feeding Tubes

MGSIL4/40	4FR x 40cm Silicone	10
MGSIL5/40	5FR x 40cm Silicone	10
MGSIL5/60	5FR x 60cm Silicone	10
MGSIL5/90	5FR x 90cm Silicone	10
MGSIL6/40	6FR x 40cm Silicone	10
MGSIL6/60	6FR x 60cm Silicone	10
MGSIL6/90	6FR x 90cm Silicone	10
MGSIL8/40	8FR x 40cm Silicone	10
MGSIL8/60	8FR x 60cm Silicone	10
MGSIL8/90	8FR x 90cm Silicone	10

Accessories

C005M	Medicine Straw	50
C020M	Milk Straw	50
MC01	Male Cap (Packed Singly)	50
MC08	Male Cap (Pack of 8)	50
MK01	Male Adapter	50
MN01M	Medicine Needle	50
MF01M	Filter Straw	50
MA01	Medicine Bottle Adapter	50
MA02	Bottle Adapter 14mm - 16mm	50
MA03	Bottle Adapter 17mm - 20mm	50
MA04	Bottle Adapter 22mm - 24mm	50
MA05	Bottle Adapter 18mm - 21mm	50

Surfactant Administration Kit

MSF01	Surfactant Kit	25
-------	----------------	----

Double Lumen Tubes

MNF02	Neonatal Replugle Tube	5
-------	------------------------	---



US Distributor:

Marian Medical, Inc.

3600 Chamberlain Lane, Suite 122, Louisville, KY 40241

Ph: 502-425-6363, email: cs@marianmedicalonline.com,

www.marianmedicalonline.com



SUPERLOCK[®]

enteral feeding system

